

RICOH
imagine. change.

Ricoh Biosciences, Inc.



Amyotrophic Lateral Sclerosis (ALS)

Amyotrophic Lateral Sclerosis (ALS)

Building Better ALS Models for Translational Discovery

ALS drug discovery has long been limited by models that fail to capture human disease biology. Our iPSC-derived ALS model recapitulates hallmark patient phenotypes such as TDP-43 mislocalization and motor neuron degeneration—without artificial overexpression. Partner with us to design, optimize, and execute custom cell-based assays that generate meaningful, translational data to drive your research forward.

TDP-43 Pathology: A Hallmark of ALS

A defining hallmark of ALS pathology is the mislocalization of TDP-43 protein from the nucleus to the cytoplasm, where it forms aggregates that disrupt RNA processing and drive motor neuron degeneration. Ricoh Biosciences' iPSC-derived ALS models reproduce this key disease feature, providing a human-relevant platform to investigate disease mechanisms and evaluate therapeutic strategies targeting TDP-43 pathology.




Figure 1: Representative immunofluorescence of iPSC-derived motor neurons from a sporadic ALS donor (bottom) and a healthy control (top). Nuclei (blue), and total TDP-43 (green).

Immunofluorescence imaging of iPSC-derived motor neurons (see Figure 4) reveals distinct TDP-43 localization patterns in ALS versus control cultures. In healthy neurons, TDP-43 remains confined to the nucleus; in ALS models, it accumulates in the cytoplasm — faithfully reproducing a defining pathological signature of ALS.

Quantifying TDP-43 in ALS Motor Neurons

Quantitative image analysis confirms a significant increase in cytoplasmic TDP-43 puncta in ALS cultures relative to healthy controls (**** $p < 0.0001$). This robust phenotype provides a reproducible assay for evaluating therapeutic compounds that modulate TDP-43 aggregation or localization dynamics in a human-relevant system.




Figure 2: Quantification of TDP-43 localization in Quick-Neuron™ motor neurons. Violin plots show the ratio of cytoplasmic to total TDP-43 puncta in cultures derived from a sporadic ALS patient and a healthy control. Significance was quantified by T-test. Stars denote statistical significance: (**** = $p < 0.0001$).

Amyotrophic Lateral Sclerosis (ALS)

Custom assay services and human, iPSC-derived ALS disease models for drug discovery.

Compound Testing in iPSC-Derived ALS Models.

To demonstrate the translational utility of our iPSC-derived ALS motor neuron models, we tested compounds targeting key disease mechanisms. Edaravone, CK1 inhibitors, and Sephin1—representing antioxidant, TDP-43 regulatory, and stress-response pathways—each produced measurable rescue of TDP-43 mislocalization in ALS motor neurons. These findings highlight the model's sensitivity to diverse therapeutic mechanisms and its value as a platform for screening candidate ALS treatments.




Figure 3: Violin plots show the distribution of per-well mean cytosolic puncta/total puncta ratios for healthy control and ALS cultures treated with vehicle (0.1% DMSO) or edaravone (10 or 100 µM). ALS motor neurons treated with DMSO show elevated cytoplasmic TDP-43, consistent with ALS pathology. 24-hour treatments with Edaravone at both 10 µM and 100 µM concentrations significantly rescues this mislocalization, reducing cytosolic TDP-43 levels.



Physiologically Relevant

Recapitulates disease hallmarks such as TDP-43 mislocalization.



Consistent & Reproducible

Reliable performance and batch-to-batch consistency.



Sporadic & Familial ALS Models

Patient-derived lines for sporadic and familial ALS.



Naturally Evolving Pathology

Progressive, spontaneous disease phenotypes.

Amyotrophic Lateral Sclerosis (ALS)




Figure 4. Quantification of the cytoplasmic-to-total TDP-43 ratio in Quick-Neuron™ motor neuron cultures derived from a sporadic ALS patient and a healthy control. All cultures were treated with DMSO, and ALS cultures were additionally treated with a CK1 inhibitor (1 μM or 10 μM).




Figure 5. Quantification of the cytoplasmic-to-total TDP-43 ratio in Quick-Neuron™ motor neurons derived from a sporadic ALS patient and a healthy control. All cultures were treated with DMSO, and ALS cultures additionally received Sephin1 (100 nM or 500 nM).

Our ALS Disease Modeling Services.



Assay Services

Reliable testing with iPSC-derived cells for drug discovery and safety evaluation.



Assay Development

Custom assay design tailored to your research goals, including disease models and endpoints.



Differentiation Services

Rapid, high-quality differentiation of iPSCs into virtually any human cell type.



Let's Connect



www.biosciences.ricoh.com



info@biosciences.ricoh.com



+1 (443) 869-5420